


Product Selector Guide



- A** Offset neck design to maintain -150°C or -190°C in vapor storage and provide low liquid nitrogen consumption
- B** Optimized lid design for enhanced durability
- C** Rotating interior tray provides easy access to cryobiological samples
- D** Low Maintenance, all-stainless steel construction
- E** Annular filling lines reduce frost and ice formation near lid
- F** Durable casters
- G** Two-level folding step on 1500 and 1800 series models provides lowest liftover height
- H** Advanced TEC 3000 controller

